

### 1.50 SEMI-FINISHED SV

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Single vision spherical	1.50	0.00, 0.50, 1.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00, 11.00, 12.00, 14.00	Clear	65	UC/HC
Single vision spherical	1.50	0.00, 0.50, 1.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00	Clear	70	UC/HC
Single vision spherical	1.50	0.00, 0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00	Clear	75	UC/HC
Single vision aspherical	1.50	0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.00	Clear	70	UC/HC
Single vision aspherical	1.50	0.00, 2.00, 3.00, 4.00, 6.00, 7.00, 8.00	Clear	75	UC/HC
Single vision spherical	1.50	2.00, 4.00, 6.00, 8.00	Polarized	75	UC/HC
Single vision spherical	1.50	2.00, 4.00, 6.00, 7.00, 8.00	Polarized	80	UC/HC
Single vision spherical	1.50	0.00, 1.00, 2.00, 3.00, 4.00, 5.25, 6.00, 7.25, 8.00, 9.00, 10.00	Clear	80	UC/HC
Single vision spherical	1.50	0.00, 2.00, 4.00, 6.00, 7.00, 8.00, 9.00	Clear	85	UC/HC

### 1.50 SEMI-FINISHED BIFOCAL

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Flat top	1.50	0.50, 2.00, 4.00, 6.00, 8.00, 10.00	Clear	70/25	UC/HC
Flat top	1.50	0.50, 2.00, 4.00, 5.00, 6.00, 7.00, 8.00, 10.00	Clear	70/28	UC/HC
Flat top	1.50	2.00, 4.00, 6.00, 8.00	Polarized	77/28	UC/HC
Flat top	1.50	0.50, 2.00, 4.00, 5.00, 6.00, 7.00, 8.00	Clear	75/28	UC/HC
Flat top	1.50	4.00, 6.00, 8.00	Clear	80/28	UC/HC
Round top	1.50	0.50, 2.00, 4.00, 6.00, 8.00, 10.00	Clear	70/24	UC/HC
Round top	1.50	2.00, 4.00, 6.00, 8.00	Clear	70/28	UC/HC
Blended	1.50	4.00, 6.00, 8.00	Clear	70/24	UC/HC
Blended	1.50	2.00, 4.00, 6.00, 8.00	Clear	70/28	UC/HC

### 1.50 SEMI-FINISHED PROGRESSIVE

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Progressive	1.50	2.00, 4.00, 6.00, 8.00	Clear	75/12	UC/HC
Progressive	1.50	2.00, 4.00, 6.00, 8.00	Clear	75/14	UC/HC
Progressive	1.50	2.00, 4.00, 6.00, 8.00	Clear	75/17	UC/HC
Progressive	1.50	2.00, 4.00, 6.00, 8.00	Polarized	75/14	UC/HC
Progressive	1.50	2.00, 4.00, 6.00, 8.00	Polarized	75/17	UC/HC