



1.60 SEMI-FINISHED SV

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Single vision spherical	1.60	0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00, 9.00, 10.25	Clear	70/75	UC
Single vision spherical	1.60	2.00, 4.00, 6.00, 8.00	Polarized	75	UC/HC
Single vision aspherical	1.60	0.50, 1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.25, 9.25, 10.25	Clear	70/75	UC
Single vision spherical	1.60	0.50, 2.00, 4.00, 5.00, 6.00, 7.00, 8.00	Gray	80	UC/HC

1.60 SEMI-FINISHED BIFOCAL

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Flat top	1.60	0.50, 2.00, 4.00, 6.00, 8.00	Clear	70/28	UC
Round top	1.60	2.00, 4.00, 6.00, 8.00	Clear	70/28	UC

1.60 SEMI-FINISHED PROGRESSIVE

TYPE	INDEX	BASE CURVE	COLOR	BLANK SIZE	COATING
Progressive	1.60	2.00, 4.00, 6.00, 8.00	Clear	75/14	UC
Progressive	1.60	2.00, 4.00, 6.00, 8.00	Clear	75/17	UC